

# MEXPRO®

ULTRA Light Laser Technology

 MEXDIODE®  
With Three Wavelengths Technology

MEXSHAPE®  
Bölgesel Vücut Şekillendirme Teknolojisi

MEXEMS®  
Elektro Mıknas Stimulator Technology

## LOOP

Thanks to the intense LED light on the head, it allows you to see even the finest hairs.

## COLD AIR HEAD

Coldest Head

With the state-of-the-art cooling system in **MexPro**, you will be able to perform your process without causing pain to your customers, even if you set the energy level of your device to the highest level.

**MEXPRO**  
ULTRA Light Laser Technology

## HYIGENIC SILICONE HEAD

**MexPro**'s heads are made to provide the best application for your customers. In this way, it does not cause any hygiene problems during

## LATEST TECHNOLOGY FROSTED HEAD

Designed by taking into account the shortcomings of other devices, our new head is designed by combining the integrated cold air system with the icing head.

**-25°**

# LIFETIME FREESHOT INSTALLATION

With the unique feature of MexPro, you now have the advantage of no longer having to pay for a shooting credit and can avoid additional costs.



## FASTEST EPILATION DEVICE

The **MexPro** cold air nozzle is designed to save you time. In this way, you will have the advantage of being able to shoot 4 shots per second



## THE QUIETEST EPILATION DEVICE

MexPro is the quietest device compared to other devices with its 55-60dB noise level during application. In this way, it does not disturb you and your customers during the application.



## THE MOST POWERFUL EPILATION DEVICE

MexPro Ultra Light Laser Technology is a highly effective hair removal technology that targets hair follicles and delivers wavelengths of 755nm - 810nm - 1064nm at 85% efficiency.



## NEW GENERATION LASER HAIR REMOVAL DEVICE

It can go up to **20**/j cm energy level to achieve more permanent results in a shorter time and you can achieve effective results even in the finest hairs with **85%** focused wavelength of **755 - 810 - 1064 nm**.

Thanks to the movable upper arm, it offers an ergonomic design by rotating in 2 axes and does not tire users. The details in its design are carefully designed for comfortable use.

### NO MORE PAINFUL EPILATORS

- › MexPro was created to provide your customers with a completely Painless hair removal experience. In this way, you will be able to get positive feedback from your customers after your sessions

### 95% MORE EFFECTIVE

- › 2023 Model burns hair without damaging the skin surface with its latest technology and precise results by burning hair follicles in one shot to get the most out of your time. In this way, you can use your time more efficiently and use more of it during the day allows you to apply the application..

### PERFECT SMOOTHNESS

- › With its state-of-the-art technology, among other devices, it is very effective even on the finest hairs and thanks to its powerful hair blasting effect, it perfectly smooths the skin surface.

**Confirmable Color Option**



Wight



Black



Pink



#### **Application in 4 Seasons**

Offers epilation regardless of the season.



#### **Ergonomic Design**

Easy to use thanks to the movable head.



#### **Risk Free Hair Removal**

Eliminate the risk of burns thanks to improved software.

### **TECHNICAL FEATURES**

Technical Features	New Generation Ultra Light Laser Flash Lamp
Light Spectrum Focus	630nm - 1200nm
Maximum Rate of Fire per Second	4 Hz
Skin Cooling Technology	Inverter Compressor Cooling System
User Display	15.6" Capacitive Ultra HD
Air Blowing Stage	10
Firing Area	7.6 cm <sup>2</sup>
Cold Air Blower Head Average Lifetime	1.000.000 Shots
Warranty Life of Shot Cold Air Blower Hood	24 Months Unlimited
Frosted Hood Average Lifespan	1.000.000 Shots
Frosted Nozzle Warranty Life	24 Months Unlimited
Maximum Joule Level	266 Joule
Maximum Loudness	60 dB
Cable and Hose Length	1,5 m
Maximum Power Consumption	5000 Watt
Electricity Requirement	220 - 230 VAC / 50-60 Hz
Warranty Period	24 Months



## NEW GENERATION DIODE LASER EPILATION DEVICE



**MexDiode** is a new generation laser hair removal device model **2023**. With a wavelength of **755 - 808 -1064nm**, this device is equipped with the original **600W** and **12 bars**, an impressive power. MexDiode, produced by Atayolu Medical & Cosmetics, works by targeting the hair follicles and instantly neutralizes the hair follicles without damaging the body surface. In this way, users can They can achieve results of up to **95%**.

- **MexDiode** has a unique feature; the ability to use **3** different wavelengths together. Thanks to this feature, it can effectively treat all types of hairs, including fine and light-colored hairs. This allows users to achieve successful results on a wide range of hairs.
- The device is equipped with advanced technologies such as ecovat and tech cooling 2 cooling system. This allows for perfect application even on individuals with darker skin tones. The cooling system ensures skin comfort while minimizing side effects. MexDiode has an impressive pulse capacity. With a total of **20,000,000** shots, it offers long-term and effective use. This shows the durability and longevity of the device.
- To summarize, MexDiode is a high performance laser hair removal device. It achieves effective results by targeting hair follicles, addresses a wide range of hairs with different wavelengths, provides flexibility for users with its suitability for various skin tones and has a large pulse capacity.

**ORIGINAL  
600W**  
and Flawless Power equipped  
with 12 Bar





**20,000,000 Shots**

It has an impressive firing capacity.



**95% Effective Result**

Provides 95% effective results by neutralizing hair follicles



**Dual Cooler System**

Ecovat and tech cooling cooling technology.



**Suitable for All Skin Types**

With 5 skin type options, it provides effective results on all skin types



**3 Wavelengths in One**

755 - 808 - 1064 nm

**TECHNICAL FEATURES**

User Display		10.4-inch Next Generation Full HD
Light Spectrum Focus		755 - 808 - 1064 nm
Number of Bars		12 Bars
Power Capacity		600 Watt (Original)
Energy Capacity		1-20 Joule/ cm <sup>2</sup>
Electricity Requirement		220-230 VAC / 50-60 Hertz
Application Area		22mm * 11mm
Min. in Sec. - Max. Pulse Rate		1 - 10 Hertz
Skin Type Option		5 Skin Type options
Bristle Type Option		4 different bristle type options
Head Cooling		Coherent Bar Sapphire Cooling System
Cooling Technology		Ekovat Air & Tech Cooling System
Warranty Period		24 Months



## NEW GENERATION MUSCLE BUILDING DEVICE



Using electromagnetic wave technology, **MexEms** intensively simulates and works the muscles for only 30 minutes, providing the same effect as if you did 25,000 sit-ups.

With **MexEms**, you can strengthen, shape and tighten your body by effectively working the muscles in the hips, abdomen, legs and arms.

- › MexEms' electromagnetic wave technology penetrates deeply into the muscles and intensively simulates the muscle fibers. In this way, the muscles develop and strengthen faster. Especially the muscles in the abdominal area are important for tightening and definition of the abdominal muscles. MexEms contributes to the appearance of a firm abdomen by effectively working the muscles in this area.
- › Similarly, the muscles in the buttocks, legs and arms are also targeted with MexEms. The glutes are important for achieving a firm and lifted look. MexEms offers an effective way to stimulate and sculpt the glutes. The muscles in the leg and arm areas can also be strengthened, shaped and tightened with MexEms.
- › However, keep in mind that each individual's body structure, musculature and goals are different. The effects of MexEms may vary from person to person.

**Strengthens,  
Sculpts and Firms  
your body**







In 30 minutes  
**25.000**  
Shuttle Effects



**4 Titles in One**

Designed for use on all parts of the body.



**Simple Operation**

Ease of design  
Trading is now easier



**Refreshed Technology**

Now more effective thanks to renewed software and hardwares

**TECHNICAL FEATURES**

Electro Magnetic Intensity	7 Tesla
Number of Headboard	4 Pieces
Power Capacity	300-3500 Watt
Maximum Power Consumption	3500 Watt
Operating Frequency	1-150 Hertz(Adjustable)
Electricity Requirement	220-230 VAC / 50-60 Hertz
Cable Length	2,5m
Warranty Period	24 Months

## NEW GENERATION Regional Slimming Device

- › **MexShape** has 6 heads. You can quickly and effectively eliminate unwanted fat, cellulite and problems such as striae on your body with its heads (**Radiofrequency, Vacuum, striae and ultrasonic**) that allow you to achieve perfect results.
- › Using a combination of ultrasonic and radiofrequency heads you get rid of **70%-80%** of fat. in this way, it allows you to completely get rid of cellulite problems and have your dream look.

### ULTRASONIC HEAD

- › The ultrasonic head allows you to burn fat in a short time by giving sound waves with reduced intensity to the area where you want the fat to melt in your body.

### VACUUM RADIO FREQUENCY HEAD

- › The vacuum feature of our vacuum radio frequency heads separates the fats and breaks their resistance by applying reverse pressure to the area. Thanks to this feature, fat tissues are reduced and loosened, accelerating their elimination from the body. Radio frequency feature tightens the skin it recovers and repairs the area. It is the application used in complaints such as cellulite.

### STRIA HEAD

- › Stria cap is applied to the area where stretch marks are located with its vacuum feature and eliminates the scars in a short time.

**Mex Shape does not  
cause any pain  
or discomfort during  
the procedure.**







After the 1st session, it provides a thinning  
between

**4-12 Cm**

**The result varies\***

\*from person to person



**6 Heads in One**

Radiofrequency, vacuum,  
striae in one



**Simple Operation**

Making transactions is now  
much easier with Ease of  
Design



**Refreshed Technology**

Now more effective thanks to  
renewed Software and  
hardware

**TECHNICAL FEATURES**

RF Vacuum Head	3 Pcs
Ultrasonic Head	2 Pcs
Stria Headboard	1 Piece
User Display	10" Capacitive HD
Radio frequency Headings	1.5 Mhz multipolar high frequency
Ultrasonic Heads	1.7 Mhz high frequency sound wave
Vacuum Power	700mm hg
Cable and hose length	1.5m
Maximum power consumption	1000W
Electricity Requirement	220 - 230 VAC / 50-60Hz
Warranty Period	24 Months



**ATAYOLU**  
MEDİKAL & KOZMETİK

**When you need us, we are there to help**

**48** hours within  
**TECHINICAL  
SUPPORT**

**ATAYOLU  
MEDİKAL**

0.322 451 01 01

Aydınlar Mahallesi  
Güneykırı Boyu Caddesi  
323A  
Seyhan /ADANA



[www.atayolu.com.tr](http://www.atayolu.com.tr)

